

# Veterinary ENT diagnostic instruments

## RELIABLE, LONG-LASTING INSTRUMENTS

The HEINE veterinary range offers instruments of the highest quality and functionality as well as HEINE XHL Xenon halogen bulbs and innovative power sources.

HEINE VET instruments are subjected to ongoing quality assurance checks. They stand up to the demands of the busiest veterinary practice and assure optimum levels of diagnosis and treatment.





HEINE Veterinary ENT diagnostic instruments



Slit Illumination Head



BETA200 F.O.  
 VET Otoscope



BETA100  
 VET Otoscope

	Housing construction		
Metal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Swivel lens – Magnification			
Glass/Magnification 1.6x	<input type="checkbox"/>		
Glass/Magnification 2.5x		<input type="checkbox"/>	<input type="checkbox"/>
Illumination			
F.O.	<input type="checkbox"/>	<input type="checkbox"/>	
Direct			<input type="checkbox"/>
Instrument size			
Professional	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Available power sources			
Battery 2.5V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rechargeable battery 3.5V	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tips			
Long, closed specula	<input type="checkbox"/>		
Slit metal specula, open	<input type="checkbox"/>		
Short, closed specula	<input type="checkbox"/>		
UniSpec VET disposable tips	<input type="checkbox"/>		
Tip adaptor	<input type="checkbox"/>		
VET reusable specula		<input type="checkbox"/>	<input type="checkbox"/>
AllSpec disposable tips (with adaptor)			<input type="checkbox"/>
Insufflation port		<input type="checkbox"/>	<input type="checkbox"/>
Page	031	033	033

\* with Adaptor

Slit Illumination Head and accessories



- ∴ **Solid metal construction.** Will last a lifetime.
- ∴ **Slit instrument head.** Simplifies the use of instruments.
- ∴ **Fiber optics illumination.** Bright, homogeneous light.
- ∴ **Multi-coated, scratch-resistant swivel lens with 1.6x magnification.** High-resolution image.
- ∴ **Simple open construction.** Easy to clean and disinfect with conventional solutions.
- ∴ **Double bayonet lock for specula.** Specula mounted firmly by double lock.

Slit Illumination Head	2.5V	3.5V
with swivel lens	<b>G-001.21.301</b>	<b>G-002.21.301</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.077</b>	<b>X-002.88.078</b>
Swivel lens	<b>G-000.21.209</b>	

Long, closed specula



For use with the Slit Illumination Head.

- ∴ **Metal construction, chrome-plated.** Easy to clean, disinfect and sterilize.
- ∴ **Double bayonet lock for specula.** Firm mounting of specula to head.

Closed specula for Slit Illumination Head	
57 mm long / 4 mm dia.	<b>G-000.21.330</b>
65 mm long / 6 mm dia.	<b>G-000.21.331</b>
90 mm long / 7 mm dia.	<b>G-000.21.332</b>
150 mm long / 10 mm dia.	<b>G-000.21.350</b>

57 mm 65 mm 90 mm 150 mm  
 4 mm dia. 6 mm dia. 7 mm dia. 10 mm dia.

Slit metal specula, open



For use with the Slit Illumination Head.

- ∴ **High-quality chrome-plated metal construction.** Easy to clean, disinfect and sterilize.
- ∴ **Matt internal surfaces.** Reduced reflexes.
- ∴ **60° opening with rounded edges.** Ideal for instrumentation, atraumatic insertion.
- ∴ **Double bayonet fitting.** Positive fit.

Slit metal specula for Slit Illumination Head	
65 mm long / 6 mm dia.	<b>G-000.21.341</b>
90 mm long / 7 mm dia.	<b>G-000.21.342</b>

65 mm long 90 mm long  
 6 mm dia. 7 mm dia.

### Short, closed specula

#### for Slit Illumination Head and BETA100 VET Otoscope



5,5 mm dia. 4,5 mm dia. 3,5 mm dia.



2,8 mm dia. 2,2 mm dia.



- ⌘ **Made of high-density SANALON S with a smooth external surface.** Dust and germs are easily removed. Disinfection and autoclaving at up to 134 °C is possible.
- ⌘ **Ergonomic shape.** Painless atraumatic insertion.
- ⌘ **Integral metal ring with 2 bayonet locks.** Firm mounting of specula to head.
- ⌘ **40 mm long.**

#### Short, closed SANALON S specula

##### for Slit Illumination Head and BETA100 VET Otoscope

1 set = 5 veterinary specula in SANALON S	G-000.21.316
40 mm long / 2.2 mm dia.	G-000.21.310
40 mm long / 2.8 mm dia.	G-000.21.311
40 mm long / 3.5 mm dia.	G-000.21.312
40 mm long / 4.5 mm dia.	G-000.21.313
40 mm long / 5.5 mm dia.	G-000.21.314

### HEINE UniSpec® disposable ear specula

#### for BETA100 and with Adaptor for Slit Illumination Head



2,5 mm dia.

4 mm dia.



[01]

The adaptor illustrated [01] is necessary when using UniSpec specula with the Slit Illumination Head.


- ⌘ **No cleaning necessary.** Saves time.
- ⌘ **Hygienic and safe.** Eliminates cross-infection.
- ⌘ **Packed in units of 1000.** Assures a competitive price.
- ⌘ **Original HEINE design and manufacture.** Guaranteed good fit to the instrument head.

#### UniSpec disposable ear specula (grey)

Pack of 1000 UniSpec disposable specula 2.5 mm dia.	B-000.11.242
Pack of 1000 UniSpec disposable specula 4 mm dia.	B-000.11.241
Adaptor for use with Slit Illumination Head [01]	G-000.21.302

⊗: Single use only.

### Spreadable speculum for Slit Illumination Head




This versatile instrument is widely used as a universal speculum for small animals.

- ⌘ **For use with the Slit Illumination Head.**
- ⌘ **Distal jaw opens up to 28 mm with 4 click-stops.** Useful for examination and therapy.
- ⌘ **High-quality, chrome-plated metal construction.** Easy to clean, disinfect and sterilize.
- ⌘ **Matt internal surfaces.** Avoids reflexes.
- ⌘ **Oval branches with rounded edges.** Ideal for instrumentation, atraumatic.
- ⌘ **2 bayonet locks.** Firm mounting of spreadable Speculum to head.

Spreadable speculum for Slit Illumination Head (95 mm long)	G-000.21.360
---	--------------

HEINE BETA®200 Fiber Optic (F.O.) VET Otoscope

**Veterinary Otoloscope with Fiber Optic illumination**



Fiber Optic (F.O.) illumination


Solid, all-metal construction. A quality instrument for the Vet's office with Fiber Optic (F.O.) illumination. Insufflation port. The swivel window features an airtight seal.

- ⋄ **Distal Fiber Optic illumination.** Perfectly even illumination of auditory canal and tympanum.
- ⋄ **Large, scratch-resistant glass viewing window with 2.5x magnification.** Unobstructed view, distortion-free image.
- ⋄ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ⋄ **Insufflation port.** Reliable test of tympanic mobility.
- ⋄ **Durable metal body, matt chrome-plated.** Unbreakable, long-life.
- ⋄ **XHL Xenon Halogen Technology.** Extremely bright, white light.
- ⋄ **For use with reusable SANALON S VET Specula (G-000.21.213).**
- ⋄ **Instrument head matt-black inside.** Eliminates reflexes.

<b>BETA200 VET Otoloscope</b>	<b>2.5V</b>	<b>3.5V</b>
without specula, without handle	<b>G-001.21.250</b>	<b>G-002.21.250</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.077</b>	<b>X-002.88.078</b>

HEINE BETA®100 VET Otoloscope

**Veterinary Otoloscope with direct illumination**



Direct Xenon Halogen illumination

Otoloscope with all-metal construction and direct Xenon Halogen illumination. Suitable for use with short, closed specula. With the use of the VET adaptor, it is possible to fit reusable VET specula for ear examinations.

- ⋄ **All-metal housing.** Exceptionally-tough and long-lasting.
- ⋄ **XHL Xenon Halogen Technology.** Extremely bright, white light.
- ⋄ **Large, scratch-resistant glass viewing window with 2.5x magnification.** Unobstructed view, distortion-free image.
- ⋄ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ⋄ **Insufflation port.** Reliable test of tympanic mobility.
- ⋄ **Can be used with short, closed specula G-000.21.316 (no adaptor required).**
- ⋄ **Versatile, useful for very small animals.**
- ⋄ **Instrument head matt-black inside.** Eliminates reflexes.

<b>BETA100 VET Otoloscope</b>	<b>2.5V</b>	<b>3.5V</b>
without specula, without handle	<b>G-001.21.260</b>	<b>G-002.21.260</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>	<b>X-002.88.049</b>
Adaptor for reusable SANALON S VET specula [01]	<b>G-000.21.214</b>	
Adaptor for AllSpec tips [02]	<b>B-000.11.306</b>	

[01]

[02]

## Reusable VET specula

### for BETA200 VET or BETA100 VET with adaptor



- ∴ **Made of high-density SANALON S with a smooth external surface.** Dust and germs are easily removed. Disinfection and autoclaving at up to 134 °C is possible.
- ∴ **Ergonomic shape.** Atraumatic.
- ∴ **Convenient fixation system.** Positive fit to otoscope head.

#### Reusable SANALON S specula

1 set = 3 reusable SANALON S specula	<b>G-000.21.213</b>
61 mm long / 4 mm dia.	<b>G-000.21.210</b>
61 mm long / 6 mm dia.	<b>G-000.21.211</b>
61 mm long / 9 mm dia.	<b>G-000.21.212</b>

## HEINE AllSpec® Disposable tips

### for BETA100 VET with adaptor



- ∴ **For single use.** Hygienic, no cross-infection.
- ∴ **Original HEINE design.** Firm fit to the otoscope head.
- ∴ **Assured quality.** No sharp edges.

#### AllSpec Disposable tips (grey)

Pack of 1000 pcs. 2.5 mm dia. [03]	<b>B-000.11.128</b>
Pack of 1000 pcs. 4 mm dia. [03]	<b>B-000.11.127</b>
Case of 10 dispenser packs of 250 tips, 2.5 mm dia [01]	<b>B-000.11.151</b>
Case of 10 dispenser packs of 250 tips, 4 mm dia [01]	<b>B-000.11.150</b>
AllSpec Tip-dispenser [02]	<b>B-000.11.148</b>



#### AllSpec Tip-dispenser [02]

**The right size, visible and instantly available.** A clean new tip for each and every patient. Contents: 50 pcs. Ø 2.5 mm and 50 pcs. Ø 4 mm.

⊗: Single use only.

## Insufflation Bulb

### For BETA200VET and BETA100VET Otoscope.



- ∴ **For pneumatic testing.**

Insufflation bulb	<b>B-000.11.240</b>
-------------------	---------------------

HEINE Veterinary Diagnostic Sets

Set complete with: **Slit illumination head**

1 set long, closed specula (57mm x 4mm dia., 65mm x 6mm dia., 90mm x 7mm dia.)  
 1 spare bulb  
 hard case



Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	<b>G-112.10.118</b>		
BETA R rechargeable handle for mains socket		<b>G-112.20.376</b>	<b>G-112.27.376</b>
BETA TR rechargeable handle with plug-in transformer		<b>G-112.20.384</b>	<b>G-112.27.384</b>
BETA NT rechargeable handle and NT300 charger		<b>G-112.20.420</b>	
BETA L rechargeable handle and NT300 charger			<b>G-112.29.420</b>

Set complete with: **Slit illumination head, BETA200 Ophthalmoscope**

1 set long, closed specula (57mm x 4mm dia., 65mm x 6mm dia., 90mm x 7mm dia.)  
 1 spare bulb for each instrument  
 hard case



Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V
BETA battery handle	<b>G-148.10.118</b>		
BETA R rechargeable handle for mains socket		<b>G-148.20.376</b>	<b>G-148.27.376</b>
BETA TR rechargeable handle with plug-in transformer		<b>G-148.20.384</b>	<b>G-148.27.384</b>
BETA NT rechargeable handle and NT300 charger		<b>G-148.20.420</b>	
BETA L rechargeable handle and NT300 charger			<b>G-148.29.420</b>

: Rechargeable handle with charge status display in the bottom insert.



## HEINE Veterinary Diagnostic Sets



Set complete with: **Slit illumination head, BETA 200 Ophthalmoscope**

5 short, closed SANALON S specula

10 UniSpec disposable ear specula with adaptor

1 set = 4 long, closed specula



(57mm x 4mm dia., 65mm x 6mm dia., 90mm x 7mm dia., 150mm x 10mm dia.)

2 slit specula

spreadable speculum

1 spare bulb for each instrument

hard case

Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V 
BETA battery handle	<b>G-150.10.118</b>		
BETA R rechargeable handle for mains socket		<b>G-150.20.376</b>	<b>G-150.27.376</b>
BETA TR rechargeable handle with plug-in transformer		<b>G-150.20.384</b>	<b>G-150.27.384</b>
BETA NT rechargeable handle and NT 300 charger		<b>G-150.20.420</b>	
 BETA L rechargeable handle and NT 300 charger			<b>G-150.29.420</b>



Set complete with: **BETA 200 Ophthalmoscope, BETA 200 VET Otoscope**

3 reusable SANALON S VET specula

1 spare bulb for each instrument

zipper case



Veterinary Diagnostic Sets	2.5V	3.5V (NiMH)	3.5V 
BETA battery handle	<b>G-151.10.118</b>		
BETA R rechargeable handle for mains socket		<b>G-151.20.376</b>	<b>G-151.27.376</b>
BETA TR rechargeable handle with plug-in transformer		<b>G-151.20.384</b>	<b>G-151.27.384</b>
BETA NT rechargeable handle and NT 300 charger		<b>G-151.20.420</b>	
 BETA L rechargeable handle and NT 300 charger			<b>G-151.29.420</b>

: Rechargeable handle with charge status display in the bottom insert.